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## D

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### Composition

The femena<sup>®</sup> / femena<sup>®</sup> gold is an intrauterine device (IUD) made of a combination of rigid polyethylene, EVA and BaSO<sub>4</sub>. The vertical stem is wrapped with copper wire. The total surface area of the copper wire is 375 mm². The femena<sup>®</sup> / femena<sup>®</sup> gold device has a width of 19.1 mm and a length of 35.5 mm. The two monofilament threads for detection and removal of the device are made of polyvinylidene fluoride. For the femena<sup>®</sup> gold version, each of the two flexible wings is additionally provided with a gold clip allowing to detect the device's position.

### Mode of action

Copper IUDs release 10-50 mg (millionths of a gram) of copper daily. This is absolutely harmless for the user since the amount of copper which is absorbed with food daily, is considerably higher. According to recent scientific research, the released copper ions interfere with energy metabolism in sperm mitochondria, which is necessary for sperm movement, both directly (competitive enzyme inhibition by Cu<sup>2+</sup> in exchange for Zn<sup>2+</sup>) and indirectly through modification of the metabolism of surface endometrial cells. Changes detected in the mucus of uterine cavities include a reduction in glucose concentration, oxygen partial pressure and pH. The migration of numerous leukocytes with their own metabolism into the mucus of the cervical canal and the uterine cavity causes a deterioration in the „biological environment“ there. The sperm swimming in the mucus are rendered immobile as they attempt to reach the ovum in the ampulla of the uterine tube. Thus, the effect of the IUD can be explained by both the prevention of fertilisation and, as a result of the action on the endometrium, the inhibition of an implantation.

### Indications

Contraception

### Instructions for use and handling

The patient shall be thoroughly informed by her physician on the risks and benefits of the femena<sup>®</sup> / femena<sup>®</sup> gold device (review of the patient package insert). She shall be advised to strictly comply with the instructions given on that insert, in particular those on warnings and self-examination.

### Contraindications

**Absolute counter-indications**

- Established or suspected pregnancy
- acute or recurrent inflammation of the upper and/or lower genital tract
- acute and recurrent inflammations of the pelvis minor or known history of recur-rent inflammatory diseases in the pelvic area
- postpartum endometritis
- uterine bleeding of unknown aetiology
- infected abortion in the previous 3 months
- uterine polyps and fibroids (myoma)
- suspected or known carcinoma of the genital tract
- untreated cervicitis and cervical dysplasia
- proven allergy to copper
- conditions with increased susceptibility to infection
- endometriosis
- positional anomalies and malformations of uterus or cervix
- history of ectopic pregnancy or if the occurrence of such pregnancy is encouraged by certain factors such as salpingitis, endometritis or peritonitis
- occurrence of a sexually transmissible disease during the last 12 months (with the exception of recurrent viral herpes, hepatitis B or cytomegalovirus infections).
- uterine length established by sounding less than 6 cm.

**Relative counter-indications:**

- secondary scars of the uterus for other reasons than a Cesarean section
- history of confirmed uterus perforation
- hyperplasia of the endometrium
- distinct dysmenorhea
- anaemia or menorrhagia
- serious heart valve defect or congenital heart diseases (The use of a IUD may in-crease the risk of a bacterial inflammation of the endocardium - endocarditis. To avoid this, the user should be treated with antibiotics when the IUD is inserted or removed.)
- long-term treatment with corticosteroids and non-steroidal antiirheumatics (see Adverse effects)
- intensive immunosuppressive therapy
- frequent change of sexual partners
- small uterine polyps and fibromas (myomas)
- short-wave therapy (pelvic therapies by means of short- and micro-waves, see Adverse effects)
- Wilson’s disease
- clotting disorders

### Warnings and precautions

Due to infection risks and their possible consequences (sterility), the insertion of the IUD in nulliparous women should be preceded by thorough evaluation of the desired advantages and the possible therapeutic risks. During and after the IUD insertion, nulliparas tend to suffer more frequently from pains than women who have already given birth. Intrauterine devices do not provide any protection neither against HIV infections nor against venereal diseases. During use of copper IUDs, heavier men-struation or dysmenorrhoea may occur. In women with very heavy menstruation, anaemia or dysmenorrhoea and women on anticoagulant medication, the femena<sup>®</sup> / femena<sup>®</sup> gold should be inserted only after carefully weighing up the advantages and disadvantages.

### Insertion of the IUD

Thorough history-taking (including vaginal swab) and an examination of the pelvic cavity are necessary in order to exclude women with contraindications. The insertion shall be made under sterile conditions by a gynecologist who disposes of the appropriate instruments. The position of the uterus and the size of the uterine cavity should be carefully determined by using a hysterosmeter in order to ensure the correct inser-tion of the IUD. The femena<sup>®</sup> / femena<sup>®</sup> gold should be inserted only where the uterine length measured by sounding is between 6 and 9 cm. The uterine length by sounding can be marked on the introducer with the aid of the cursor. The femena<sup>®</sup> / femena<sup>®</sup> gold is then carefully inserted until the cursor comes into contact with the cervix. The tip of the IUD has now reached the fundus of the uterus. The applicator is then with-drawn. The threads are trimmed to an intravaginal thread length of 1–2 cm. The initial follow-up examination usually takes place 4 weeks after insertion or after the first menstrual period. Subsequent follow-up examinations are to be performed every six months. Inspections include in particular the positioning check by ultrasonography. It is advisable to insert the IUD during or shortly after menstruation. There is then no need for dilatation. Insertion of the IUD is also possible immediately after an artificial abortion in the first trimester. An IUD may be inserted at the earliest 6 weeks postpar-tum. The risk of perforation of the uterus and expulsion is increased after an abortion or a birth. After a caesarean section, it is recommended that 12 weeks should elapse before insertion of an IUD. For insertion of the IUD, see figures 1 to 3

### Removal of the IUD

The femena<sup>®</sup> / femena<sup>®</sup> gold should be removed at the latest after 5 years, prefe-rably during menstruation, and replaced as necessary with a new device. The two threads on the femena<sup>®</sup> / femena<sup>®</sup> gold grasped with forceps as close as possible to the external of uteri and the IUD is removed by carefully pulling on the threads. After removal, the femena<sup>®</sup> / femena<sup>®</sup> gold must be checked for completeness. Any retained fragments may possibly be expelled painlessly during menstruation. If at-tempted removal of an IUD is unsuccessful or if a fragment is retained, diagnostic investigations should be performed to establish whether perforation or embed-ment has occurred. The IUD contains small amounts of the contrast medium barium sulphate and is therefore radiopaque. The gold clips of femena<sup>®</sup> gold version allow to run the sonographic positioning check more easily and to detect any possible fragments. Hysteroscopy also offers an opportunity for removal of fragments. Curet-tage may also be performed for removal of the IUD or of fragments. Slight pain and bleeding may occur on removal of the IUD. Precautions in case of immediate removal (see Undesirable effects).

### Undesirable effects

Basically, the following possible side effects are reported for using a copper IUD. When the IUD is inserted or immediately after such insertion, some mild cramp-ing or, in more rare cases, pains or slight dizziness may occur temporarily. However, these phenomena fade away very quickly as experience shows. For women suffering from epilepsy, insertion of the IUD may cause a sudden seizure. Accordingly, the user should be subject to a short-term surveillance treatment.

– **Pelvic infections / bleeding**

The incidence of pelvic infections in users of copper IUDs is highest in the first month after insertion and declines thereafter. Multiple sexual partners and low age are known risk factors for pelvic inflammation. As a result of a pelvic infection, fertili-ty may be impaired and the risk of an ectopic pregnancy may be increased. Pelvic inflammatory disease experienced with an IUD in situ should be treated without delay. Users must therefore be told to consult their doctor immediately if signs or symptoms are suspected. If the user experiences recurrent endometritis or repeated pelvic infections and acute infections that fail to respond to treatment within a few days, the IUD must be removed. Microbiological tests are indicated even for minor symptoms indicative of inflammatory changes. It is advisable to keep the patient under observation. The monthly menstrual period may be longer than usual. In-termediate bleeding may occur or monthly bleeding may be heavier. In some rare cases, this bleeding trouble may lead to an anaemia. In case of excessive bleeding or persistent abdominal pain of the spasmodic type, it should be determined whether this is caused by the IUD. If so, the IUD must be removed.

– **Expulsion**

The first check of the position of the IUD takes place four weeks after its insertion or shortly after the first menstrual period post insertion by means of ultrasound. Partial or complete expulsion of the IUD may be accompanied by bleeding and pain. If the detection filaments are no longer palpable without an expulsion having been noticed, an examination should be performed to exclude any perforation. Ultrasonic or radiographic methods allow to localize the IUD. If an IUD is partially expelled, its contraceptive efficacy is reduced. The user should palpate the cervix at the end of each menstrual period in order to make sure that no solid plastic elements protrude into the vagina (dislocation). An IUD that has become displaced must be removed and replaced.

– **Perforation**

In case of a suspected perforation of the uterine body or cervix during insertion, the IUD must be removed immediately. In case of such a perforation on a woman who is already wearing the IUD, the device should be removed as soon as possible. Where there is any doubt, appropriate diagnostic measures should be instituted. For exclusively breastfeeding women who have a copper IUD fitted for contraception, the perforation risk is slightly higher.

– **Pregnancy**

Should a patient become pregnant with an IUD in situ, it is essential to exclude the possibility of an ectopic pregnancy. Should the woman wish to carry the pregnancy to term, careful monitoring is essential. The IUD must be removed before the end of the first trimester. The woman should be provided with detailed information on potential risks. There is no evidence to date that an IUD has caused congenital abnormalities.

### Adverse effects

On the basis of experience to date, interference with the contraceptive action of an IUD as a result of drug interactions is highly unlikely. It is reported that some medi-cations such as immunosuppressives, corticosteriods or non-steroid antiirheumatics including painkillers, may have an adverse effect on copper IUDs. If a long-term treatment with these medications is required, another contraceptive method may be more convenient. There are also indications that tetracycline antibiotics may affect proper functioning of an IUD. Therapies using high-frequency pelvic radiation should

be applied carefully on women who have a IUD fitted, since it is assumed that the copper might be subject to considerable heating inside the uterus, thus creating the danger of burns. It is reported that diathermic treatments should be considered to be safe if the applied rate does not go beyond the therapeutic limit. Due to the energetic condition of copper and according to today’s knowledge, the IUD is not subjected to the influence of nuclear magnetic resonance (NMR), magnetic resonance imaging (MRI), radiographic diagnosis or balneological therapies.

## PHARMACEUTICAL PARTICULARS

<b>List of excipients</b>	
femena <sup>®</sup>	femena <sup>®</sup> gold
Copper	Copper
Polyethylene	Gold
EVA	Polyethylene
Barium sulphate	EVA
	Barium sulphate

### Incompatibilities

No known incompatibilities.

### Shelf-life

5 years

### Special precautions for storage

Store at 15-30 °C. Protect from direct sunlight and moisture.

### Nature and contents of container

Intrauterine devices and accessories are supplied in a sealed sterile envelope made of composite film to DIN 868.

### Instructions for use and handling

The femena<sup>®</sup> / femena<sup>®</sup> gold is individually packaged and sterilised (gas sterilisa-tion). For single use only. Not to be resterilised. Not to be used if the packaging is damaged or opened. Not to be used after the expiry date printed on the label of the packaging. IUDs are always to be inserted under aseptc conditions. On removal, the femena<sup>®</sup> / femena<sup>®</sup> gold is to be disposed of in accordance with the latest guidelines on handling and disposal of potentially infectious materials.

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HA 2 Medizintechnik GmbH
HEXAL AG