Hyaluron HEXAL®

DE

CZ

GB

SK

DK

FI

NO

NL

SE

PL

Hyaluron HEXAL®

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Hersteller:

Lifecore Biomedical, LLC
3515 Lyman Boulevard
USA-Chaska, MN 55318

Bevollmächtigter:

HEXAL AG

Industriestraße 25
D-83607 Holzkirchen
Telefon: +49 8024 908 0
Telefax: +49 8024 908 1290
E-Mail: medwiss@hexal.com

Hyaluron HEXAL®

Sodium hyaluronate 1 %

Viscoelastic solution for injection into the joint cavity

Hyaluron HEXAL is a single use, sterile, non-pyrogenic, 1 % solution of high molecular weight sodium hyaluronate. The sodium hyaluronate in Hyaluron HEXAL is a highly purified extract of a bacterial (*Streptococcus pyogenes*) fermentation.

Composition

1 ml isotonic solution for injection contains 10 mg sodium hyaluronate, sodium chloride, monobasic sodium phosphate, dibasic sodium phosphate, and water for injection.

Indications

Hyaluron HEXAL is indicated for pain and restricted mobility due to degenerative and traumatic changes in synovial joints such as the knee.

Contraindications

Do not administer Hyaluron HEXAL to patients with known hypersensitivity (allergy) to sodium hyaluronate, hyaluronan preparations or to one of the other constituents.

Do not administer Hyaluron HEXAL to patients with cases of infections or skin diseases in the area of the injection site.

Do not administer Hyaluron HEXAL to patients if the synovial joint is severely inflamed or infected.

Do not administer Hyaluron HEXAL where joint replacement is not feasible or appropriate.

Do not administer Hyaluron HEXAL into a joint that has undergone endoprosthetic replacement.

Precautions

Caution should be exercised in patients with known hypersensitivity to drugs.

General precautions for intra-articular injections should be observed. As with any invasive joint procedure there is the rare possibility of infection.

Do not concomitantly use injection site disinfectants containing quaternary ammonium salts.

Do not inject Hyaluron HEXAL intravascularly as it may cause adverse systemic events.

Do not inject Hyaluron HEXAL extra-articularly or peri-articularly.

Hyaluron HEXAL should be used with caution in patients with evidence of lymphatic or venous stasis in the leg.

The effectiveness of a single treatment cycle of less than 3 injections of Hyaluron HEXAL has not been established.

The safety and performance of Hyaluron HEXAL when used concomitantly with other drugs or devices has not been established.

Do not use Hyaluron HEXAL if the prefilled syringe or sterile package is opened or damaged. Store in the original packaging at room temperature 2-25 °C.

Do not freeze!

If product is stored refrigerated, equilibrate product to room temperature prior to use.

zww.Verfugung.gert.Hilfres.Gebrauchs.info
If the product is re-used, the functionality and sterility of the product cannot be guaranteed. Re-use could lead to patient injury and/or infection. This product is labeled as "single use" which indicates the device is intended to be used only once for a single patient.

Do not use after the expiry date indicated on the box and labels.

Keep out of the reach of children!

Use in specific patient populations

No clinical evidence is available on the use of Hyaluron HEXAL in patients with inflammatory joint diseases such as rheumatoid arthritis or Bechterew's disease. Treatment with Hyaluron HEXAL is not recommended in these cases.

Patients with autoimmune or connective tissue disorders may be more predisposed to adverse reactions than the general population.

Pediatrics: The safety and effectiveness of Hyaluron HEXAL have not been demonstrated in children.

Pregnancy: The safety and effectiveness of Hyaluron HEXAL have not been established in pregnant women.

Lactation: It is not known whether Hyaluron HEXAL is excreted in human milk. The safety and effectiveness of Hyaluron HEXAL have not been established in lactating women.

Side effects

During or after the injection of Hyaluron HEXAL, side effects including, but not limited to pain, flushing, reddening and swelling - possibly with effusion requiring puncture (clear to clouded yellowish) - may occasionally occur in the area of the treated joint. These incidents more commonly occur in patients with already existing joint effusion and may necessitate treatment to be discontinued.

In rare cases, local and general hypersensitivity reactions such as fever, chills, anxiety, oedema, skin reactions, dyspnoea, palpitation and drop in blood pressure - in isolated cases also severe anaphylactic reactions - are possible. In these cases Hyaluron HEXAL is to be discontinued immediately; the necessary symptomatic emergency measures are to be initiated.

Interactions with other products

No information on the compatibility and effectiveness of Hyaluron HEXAL with other solutions for intra-articular use is available to date. The concomitant use of an oral analgesic or anti-inflammatory drug during the first few days of treatment may be helpful for the patient.

Dosage and administration

Hyaluron HEXAL should only be administered by a physician. Strict aseptic administration technique must be followed. Inject Hyaluron HEXAL into the affected joint once a week for a total of 3-5 injections. Several joints may be treated at the same time. Depending on the severity of the joint disease, the beneficial effects of a treatment cycle of three to five intra-articular injections may last up to six months. Treatment cycle may be repeated upon necessity. In case of joint effusion it is advisable to reduce the effusion prior to injection. Treatment with Hyaluron HEXAL can be started two or three days later. As with any invasive joint procedure, it is recommended that the patient avoids any strenuous activities or prolonged (i.e., more than 1 hour) weight-bearing activities such as jogging or tennis within 48 hours following the intra-articular injection.

The contents and outer surface of the Hyaluron HEXAL pre-filled syringe are sterile as long as the pouch remains unbroken. Take the pre-filled syringe out of the pouch, remove the Luer Lok cap from the syringe, attach a suitable needle (for example 19 or 21 G) and secure it by turning slightly. Remove the air bubble from the syringe prior to injection. Inject the full volume of Hyaluron HEXAL into the treated joint.

Characteristics and mode of action

Hyaluronan is found throughout the tissues of the body with high concentrations in the synovial fluid, vitreous humor and skin. Synovial fluid, which is viscoelastic due to the presence of hyaluronan, is found in all synovial joints, particularly the large weight-bearing joints, where it ensures normal, painless movement due to its lubricating and shock-absorbing properties.

In degenerative joint disorders such as osteoarthritis, the viscoelasticity of the synovial fluid is markedly reduced thereby decreasing its lubricating and shock-absorbing functions. Supplementing this synovial fluid with intra-articular injections of highly purified hyaluronan can ameliorate some of the symptoms of osteoarthritis. A decrease in pain and an improvement in joint mobility may last for several months after a treatment cycle of three to five intra-articular injections.

Instructions for disposal

After application, syringe, needle and any material not used are to be discarded in accordance with the regulations of the country concerned.

Presentation and package size

One pre-filled syringe of 20 mg/2 ml Hyaluron HEXAL in a sterile pack. Filter sterilized. Secondarily sterilized by ethylene oxide.

To be used only by a physician

Manufactured by:
Lifecore Biomedical, LLC
3515 Lyman Boulevard
USA-Chaska, MN 55318

Authorized Representative: **HEXAL AG**

Industriestraße 25
D-83607 Holzkirchen
Phone: +49 8024 908 0
Telefax: +49 8024 908 1290
E-mail: medwiss@hexal.com

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