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## Handgelenk **MOBIL SLIM**

**PAD**  
ARRHYTHMIE  
ERKENNUNG

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BENUTZER

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DE

EN

FR


RU

TR

technology by  
**microlife<sup>®</sup>**

## aponorm® Mobil Slim – Kurzanleitung

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- ① On/Off switch
- ② Display
- ③ Battery Compartment
- ④ Cuff
- ⑤ M-Button (Memory)
- ⑥ Time Button
- ⑦ User Selection
- ⑧ LED Positioning Indicator

## Display

- |                              |                          |
|------------------------------|--------------------------|
| ⑨ Systolic Value             | ⑰ Pulse Indicator        |
| ⑩ Diastolic Value            | ⑱ Battery Display        |
| ⑪ Pulse Rate                 | ⑲ Arm Movement Indicator |
| ⑫ Date/Time                  | ⑳ Cuff Check Indicator   |
| ⑬ User Indicator             |                          |
| ⑭ Traffic Light Display      |                          |
| ⑮ Stored Value               |                          |
| ⑯ Heart Arrhythmia Indicator |                          |

Dear Customer,

Your new **aponorm**® by microlife blood pressure monitor is a reliable medical instrument for taking measurements on the wrist. It is simple to use, accurate and is highly recommended for blood pressure monitoring at home. This monitor was developed in collaboration with physicians, and clinical tests prove its measurement accuracy to be very high\*.

Please read these instructions carefully to make sure you understand all functions and safety information. We want you to be happy with your **aponorm**® by microlife product. If you have any questions, problems or want to order spare parts please contact the **aponorm**® by microlife Customer Service. Your pharmacy will be able to give you the address of the **aponorm**® by microlife distributor in your country. Alternatively, visit the internet at [www.aponorm.de](http://www.aponorm.de) where you will find a wealth of invaluable information on our products.

Stay healthy! **aponorm**® by microlife.

\* This monitor was tested in accordance with the protocol of the European Society of Hypertension (ESH) and in accordance with ISO 81060-2:2013.

Important! Please observe the information in this booklet in order to prevent damages to the blood pressure monitor and avoid error messages.

Type BF applied part

Protect from moisture

Read the instructions carefully before using this monitor.

## Table of Contents

<b>1. Important Facts about Blood Pressure and Self-Measurement</b>	20	<b>8. Battery Indicator and Battery Change</b>	
How do I evaluate my blood pressure?	21	Batteries almost flat	26
<b>2. Using the Monitor for the First Time</b>		Batteries flat - replacement	26
Inserting the batteries	22	Which batteries and which procedure to use?	27
Setting the date and time	22	Using rechargeable batteries	27
<b>3. Taking a Blood Pressure Measurement using this Monitor</b>	22	<b>9. Error messages and Troubleshooting</b>	28
Checklist for taking a reliable measurement	22	<b>10. Safety, Care, Accuracy Test and Disposal</b>	
<b>4. Display of the Heart Arrhythmia Indicator for early Detection</b>	24	Safety and protection	29
<b>5. Selection of User</b>	25	Cleaning the monitor	30
<b>6. Traffic Light Indicator in the Display</b>	25	Cleaning the cuff	30
<b>7. Data Memory</b>		Accuracy test	30
Viewing the stored values	25	Disposal	30
Memory full	25	<b>11. Warranty</b>	31
Clearing all values of an individual user	26	<b>12. Technical data</b>	32
Not storing a reading	26	Warranty card (see back cover)	

## Important Facts about Blood Pressure and Self-Measurement

- **Blood pressure** is the pressure of the blood flowing in the arteries generated by the pumping of the heart. Two values, the **systolic** (upper) value and the **diastolic** (lower) value, are always measured.
- The monitor indicates the **pulse rate** (the number of times the heart beats in a minute).
- **Permanently high blood pressure values can damage your health and must be treated by your doctor!**
- Always discuss your values with your doctor and tell them if you have noticed anything unusual or feel unsure. **Never rely on blood pressure readings alone.**
- **A number of factors can affect the accuracy of measurements taken on the wrist.** In some cases the results may differ from the measurements taken on the upper arm. We therefore advise you to compare these values with those produced by the upper arm measurement and discuss them with your doctor.
- There are many causes of **excessively high blood pressure values.** Your doctor will explain them in more detail and offer treatment where appropriate. Besides medication, weight loss and exercise can also lower your blood pressure.
- **You must not alter the dosages of any drugs prescribed by your doctor under any circumstances!**
- Depending on physical exertion and condition, blood pressure is subject to wide variations as the day progresses. **You should therefore take your measurements in the same quiet conditions and when you feel relaxed!** Take at least two measurements per day, at the same time every day, one in the morning and one in the evening. Take your medication **after** the measurement.
- It is quite normal for two measurements taken in quick succession to produce **significantly different results.**
- Even considerable **deviations** between measurements taken by your doctor or in the pharmacy and those taken at home are quite normal, as these situations are completely different.
- **Several measurements** provide a much clearer picture than just one single measurement.
- Leave a small break of at least 15 seconds between two measurements.

- If you suffer from an **irregular heartbeat** (arrhythmia, see „Section 4.“), measurements taken with this monitor should be evaluated by your doctor.
- **The pulse display is not suitable for checking the frequency of pacemakers!**
- If you are **pregnant**, you should monitor your blood pressure very closely as it can change drastically during this time!

### How do I evaluate my blood pressure?

Table for classifying home blood pressure values in adults in accordance with the international Guidelines (ESH, AHA, JSH).  
Data in mmHg.

	Range	Systolic	Diastolic	Recommendation
	Blood pressure too low	<100	<60	Consult your doctor
1.	Blood pressure optimum	100 - 130	60 - 80	Self-check
2.	Blood pressure elevated	130 - 135	80 - 85	Self-check
3.	Blood pressure too high	135 - 160	85 - 100	Seek medical advice
4.	Blood pressure far too high	>160	>110	Urgently seek medical advice

The higher value is the one that determines the evaluation. Example: a blood pressure value of 140/80 mmHg or a value of 130/90 mmHg **indicates** „blood pressure too high“.

## 2. Using the Monitor for the First Time

### Inserting the batteries

After having unpacked the monitor, insert the batteries. The battery compartment ③ is on the back of the monitor. Insert the batteries (2 x size AAA 1.5 V), observing the indicated polarity.

### Setting the date and time

1. After the new batteries are fitted, the year flashes in the display. Press the M button ⑤ to set the year. Press the time button ⑥ for confirmation and for changing to the month setting mode.
2. Now the correct month can be set by pressing the M button. Press the time button for confirmation and for changing to the date setting mode.
3. Follow the instructions above to set the date, hour and minutes.
4. Once you have set the minutes and pressed the time button ⑥, the date and time are set, and the time is displayed.
5. If you want to change the date and time again, press and hold the time button down for approx. 3 seconds until the year starts to flash. Now you can enter the new values as described above.

## 3. Taking a Blood Pressure Measurement using this Monitor

### Checklist for taking a reliable measurement

1. Avoid activity, eating or smoking immediately before taking the measurement.
2. Sit down for at least 5 minutes before the measurement and relax.
3. Always take the measurement when seated and always on the same wrist. Take your measurements on the wrist where the blood pressure is normally higher. Do not cross your legs.
4. Remove any clothing on your wrist and your watch so that your wrist is free.
5. Always ensure the cuff is positioned correctly, as shown in the pictures in the **short instruction card**.
6. Fit the cuff closely, but not too tightly. The cuff will cover a wrist diameter of 13.5 to 23 cm (5.25 - 9,05 inches). The monitor comes to sit on the inside of your wrist and must be at the same height as your heart.
7. **Select a user by setting the user switch ⑦ to either 1 or 2.**

8. Press the ON/OFF button ① to start the measurement.
  9. The LED positioning indicator ⑧ switches on automatically. Support your arm in a relaxed position at about the same height as your heart and slightly move your arm until the LED light is brightest and located in the centre of the positioning indicator.
  10. The measurement will begin after 3 seconds, and the cuff will inflate automatically. It is completely normal that the pumping speed slows down. Relax, do not move and do not tense your arm muscles until the measurement result is displayed. Breathe normally and do not talk.
  11. The measurement is performed while the cuff inflates. The pumping speed can vary, which is completely normal.
  12. During the measurement the pulse indicator ⑰ in the display is flashing.
  13. The result, comprising the systolic ⑨ and the diastolic ⑩ blood pressure and the pulse rate ⑪ is displayed. Note also the explanations on further display symbols in this booklet.
  14. Remove and switch off the monitor (the monitor will switch off automatically after approx. 1 minute) and enter the result in the enclosed blood pressure pass.
-  You can stop the measurement at any time by pressing the ON/OFF button (e.g. if you feel unwell or the pressure of the cuff becomes unpleasant).

## 4. Display of the Heart Arrhythmia Indicator for early Detection

This symbol ⑩ indicates that certain pulse irregularities were detected during the measurement. If the result deviates from your normal blood pressure repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

### **Information for the physician in the event of frequent display of the arrhythmia indicator**

This device is an oscillometric blood pressure monitor that also analyses the pulse frequency during measurement.

The monitor was subjected to a clinical study.

The arrhythmia symbol is displayed after a measurement, if pulse irregularities have occurred during this measurement.

If the symbol appears frequently (e.g. several times per week with measurements performed daily) we recommend the patient to seek an in-depth diagnostic work-up.

The use of this monitor does not replace a cardiac examination, but helps detect pulse irregularities at an early stage.

## 5. Selection of User

This monitor stores the data of 2 different users.

Before each measurement slide the user switch ⑦ in the position of the user in question.

User 1: Slide the user switch ⑦ to the left to the „User 1“ icon.

User 2: Slide the user switch ⑦ to the right to the „User 2“ icon.

☞ The person performing the first measurement should select „User 1“.

## 6. Traffic Light Indicator in the Display

The bars on the left-hand edge of the display ⑭ show the range within which the indicated blood pressure value lies. Depending on the height of the bar, the readout value is either within the optimum (green), elevated (yellow), too high (orange) or dangerously high (red) ranges. The classification is based on the 4 ranges in the table as defined by the international guidelines (ESH, AHA, JSH), as described in „Section 1.“.

## 7. Data Memory

At the end of a measurement, this monitor automatically stores each result, including date and time.

### Viewing the stored values

Select user 1 or 2 by shifting the user switch ⑦. Press the M button ⑤ briefly, when the monitor is switched off. The display first shows „M“ and a figure, e.g. „M 17“. This means 17 values are stored in the memory. The monitor then switches over to the last result stored. Pressing the M button again displays the previous value. Pressing the M button repeatedly allows you to move from one stored value to another.

### Memory full

☞ Please note that the maximum memory capacity is 60 values. When the 60 value memory is full, the oldest value is automatically overwritten by the 61st value. Values should be evaluated by a doctor before the memory capacity is reached and the data are overwritten and thus lost.

## 8. Battery Indicator and Battery Change

### Clearing all values of an individual user

Ensure that you selected the intended user. If you are sure that you want to permanently remove all stored values, hold down the M button (the monitor must have been switched off beforehand) until „CL“ appears and then release the button. To permanently clear the memory, press the M button while „CL“ is flashing. It is not possible to clear selected values.

### Not storing a reading

As soon as the result appears in the display, press the ON/OFF button ① and keep it pressed until „M“ ⑮ is flashing. Confirm the deletion by pressing the M button ⑤ again.

### Batteries almost flat

When the batteries are exhausted to approximately 3/4 of their full charge, the battery symbol ⑱ will flash as soon as the monitor is switched on (only part of the battery charge left). Although the monitor will continue to measure reliably, you should obtain replacement batteries.

### Batteries flat – replacement

When the batteries are flat, the flat battery symbol ⑲ will flash as soon as the monitor is switched on. You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment ③ by sliding the cover away.
2. Replace the batteries – ensure correct polarity as shown by the symbols in the compartment.
3. To set date and time, follow the procedure described in „Section 2“.

☞ values are retained in the memory although date and time must be reset; therefore the year flashes automatically after battery replacement.

### **Which batteries and which procedure to use?**

- ☞ Please use 2 new, long-life 1.5 V, size AAA alkaline batteries. Use only batteries of one brand at a time.
- ☞ Do not use batteries beyond their date of expiry.
- ☞ If the monitor is not going to be used for a prolonged period of time remove the batteries.

### **Using rechargeable batteries**

You can also operate this monitor with rechargeable batteries.

- ☞ Please use only „NiMH“ type rechargeable batteries.
- ☞ When the flat battery symbol appears, remove the batteries and recharge them. They must not remain in the monitor as they may become destroyed (exhaustive discharge due to low consumption of the monitor even when switched off).
- ☞ Remove the rechargeable batteries from the monitor if you do not intend to use it for one week or more.
- ☞ The batteries canNOT be charged while inserted in the blood pressure monitor. Recharge the batteries in an external charger and observe the information regarding charging, handling and service life.

## 9. Error Messages and Troubleshooting

When an error occurs during a measurement, the measurement is interrupted and an error message displayed, e.g. „ERR 3“.

Error	Designation	Potential cause and remedy
«ERR 1»	Signal too weak	The pulse signals on the cuff are too weak. Re-position the cuff and repeat the measurement.
«ERR 2» ①9	Error signal	During the measurement interference signals were detected on the cuff due to, for instance movement or muscle tension. Repeat the measurement, not moving your arm.
«ERR 3» ②0	No pressure in the cuff	It is not possible to build up sufficient pressure in the cuff. There may be a flaw in the cuff, causing a pressure leak. Check that the cuff connection is correct and that the cuff is not too loose around your wrist. Replace the batteries if necessary. Repeat the measurement.

Error	Designation	Potential cause and remedy
«ERR 5»	Abnormal result	The measuring signals are inaccurate and no result can therefore be displayed. Read through the checklist for performing reliable measurements and repeat the measurement.
«HI»	Pulse or cuff pressure too high	The pressure in the cuff is too high (over 300 mmHg) or your pulse rate is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
«LO»	Pulse too low	The pulse rate is too low (less than 40 beats per minute). Repeat the measurement.*

\* Please consult your doctor, if this or any other problem occurs repeatedly.

☞ If you consider the measurement results to be unusual, please observe the information in „Section 1.“.

## 10. Safety, Care, Accuracy Test and Disposal

### Safety and protection

- This monitor may not be used for purposes other than those described in these instructions. The manufacturer will not accept any liability for damage caused by improper use.
- This blood pressure monitor contains sensitive components and must be handled carefully. Observe the storage and operating conditions described in the „Technical data“ section.
- Protect the monitor from:
  - water and moisture
  - impact and falls
  - contamination and dust
  - direct sunlight
  - heat and cold
- The cuffs are rather delicate and must be handled carefully.
- Inflate the cuff only when fitted to your wrist.
- The function of this monitor can be impaired by strong electromagnetic fields as emitted by mobile telephones or radio installations. We recommend observing a minimum distance of 1 metre. If you cannot maintain the minimum distance please check the proper function of the monitor before use.

- Do not use the monitor if you think it is damaged or if you notice anything unusual.
- Never open the monitor.
- If the monitor is not going to be used for a prolonged period of time the batteries must be removed.
- Please observe the additional safety information in the different sections of these instructions.

Do not allow unsupervised children to use the monitor. Some of the components are so small that children can swallow them or choke on them.

### **Cleaning the monitor**

Clean the monitor with a soft, dry cloth only.

### **Cleaning the cuff**

Stains on the cuff can be removed by careful rubbing with a damp cloth and soapsuds.

### **Accuracy test**

We recommend having this monitor tested for accuracy every 2 years or after mechanical impact (e.g. the device having been dropped). Please contact your distributor (see foreword).

### **Disposal**

Batteries and electronic devices must not be disposed of as normal domestic trash. They must be disposed of in accordance with the applicable local provisions.

## 11. Warranty

This monitor comes with a 5 year warranty from the date of purchase. The warranty is valid only on presentation of the warranty card completed by the dealer (see back) confirming date of purchase or the receipt.

- Opening of or modification to the monitor by the user voids all warranty claims.
- The warranty does not cover damage caused by improper handling, battery leakage, accidents or failure to observe the instructions for use.
- The warranty covers the monitor and its cuff. Battery and wear parts are not covered by the warranty.

Please contact your distributor (see foreword).

## 12. Technical data

### Operating conditions:

10 - 40 °C / 50 - 104 °F

15% to 95% relative maximum humidity of air

### Storage conditions:

(-20) – (+55) °C / (-4) – (+131) °F

15% to 95% relative maximum humidity of air

**Weight:** 115 g

**Dimensions:** 81 x 65 x 21 mm

### Measuring method:

oscillometric, corresponding to the Korotkoff method:  
phase I systolic, phase V diastolic

### Measuring range:

20 - 280 mmHg – blood pressure

40 - 200 beats per minute – pulse

**Display range of cuff pressure:** 0 – 299 mmHg

**Measuring resolution:** 1 mmHg

**Statistic accuracy:** pressure within  $\pm 3$  mmHg

**Pulse accuracy:**  $\pm 5$  % of the readout value

**Voltage source:** 2 x 1,5 V Alkaline-Batterien, Größe AAA

**Battery life:** approx. 320 measurements (new batteries)

**IP class:** IP20

### Reference to standards:

EN 1060-1 /-3 /-4; IEC 60601-1;

IEC 60601-1-2 (EMC); IEC 60601-1-11

### Average useful life:

Monitor: 5 years or 10,000 measurements

Accessories: 2 years

This monitor complies with the requirements of the Medical Device Directive 93/42/EEC.

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Microlife AG
WEPA Apothekenbedarf GmbH & Co KG