# MEDISANA®

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## 1.1 Thank you very much

Thank you for your confidence in us and congratulations on your purchase!  
The electric lice comb **LCS** you have purchased is a high-quality product manufactured by **MEDISANA**.

To achieve the best results and ensure a long life of reliable operation of your **MEDISANA** electric lice comb **LCS**, we recommend that you carefully read through the following instructions on using and caring for the product.

## 1.2 Important

**Read the following safety information carefully before using the device and keep these instructions in a safe place for later use.**

- Use the unit only as intended in accordance with the directions.
- Using the unit for any other purpose will void your guarantee.
- The unit was not designed for commercial use.
- Do not leave the unit unattended near children. electric devices are not toys.
- The unit is not suited for children under the age of 3.
- If children are to use the electric lice comb, they must be supervised by an adult.
- If you have a pacemaker, artificial limbs or an electronic implant, please consult your physician before using the unit.
- The unit is not recommended for persons who suffer from epilepsy or heart disease.
- Use the unit only on dry and combed hair.
- Do not use the unit on dry, cracked or wounded skin.
- **When using the electric lice comb, avoid direct contact with your ears, eyes and mouth. Touching the metal teeth can inflict a painful, yet harmless, electrical charge.**
- The comb has rounded pearl-shaped tips on every second tooth **4**. Make sure that no teeth are missing or damaged before you use the unit. Otherwise the lice comb is no longer useable.
- In the event of a malfunction, do not attempt to repair the unit on your own; this will void your guarantee. Have all repairs carried out by an authorised service centre.
- You are only permitted to clean the unit.
- Keep the unit away from hot surfaces and guard it against exposure to direct sunlight.
- Protect the unit from moisture. Should liquid get into the unit, remove the batteries immediately and refrain from using the unit. In such a case, contact your dealer or notify us directly. Our address and phone number are located on the address page.
- Do not use the unit if the housing or comb teeth are damaged, if the unit does not function properly or if it has fallen down or fallen into water.
- Prevent accidents or damage to the unit by keeping the unit away from clothing, bristles, wires, strings, shoelaces, threads etc. when it is switched on.

## 1.2 Important

- The unit may occasionally leave your scalp tingling. This is perfectly normal and is in no way dangerous.
- Clean your lice comb after every use, in particular when more than one person is being treated, to avoid the risk of transmitting lice.
- Always switch off the unit when not in use.
- When not in use, always put the protective cap **5** back on the unit.

## 2.1 Items supplied and packaging

Please ensure that the package contents are complete. These are as follows:

- 1 **MEDISANA** electric lice comb **LCS** with cleaning brush
- 1 battery 1.5 V (type AA, LR 6)
- 1 Operating instructions

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

### WARNING

**To avoid risk of suffocation, keep packaging film away from children!**

## 2.2 What are head lice?

Head lice are small, wingless insects that live off blood on human scalps. Fully adult lice are typically 2 - 3 mm in length. Their colour varies from off-white to brown. The lice eggs are beige-brown in colour and stick to the hairline on your scalp. Head lice hatch after 7 to 10 days; the white shells they leave behind are known as nits.

An average affliction is about 10 lice. If left untreated, the number of lice on your scalp can quickly increase to 200. Lice can neither jump nor fly; this is why they can only be transmitted through direct head-to-head contact. However, it is also possible for lice to survive on textiles for a short period of time. Anyone can get lice, regardless of age, race and social standing. Lice are a problem that effect society at large; they are not just limited to schools and kindergartens. Contrary to popular belief, head lice can inhabit all types of hair – not just unwashed or dirty hair.

## 2.3 Special features of the **MEDISANA** electric lice comb **LCS**

The **MEDISANA** electric lice comb **LCS** effectively eliminates head lice without having to resort to chemical means. It works by creating an electric charge between the teeth that kills off the lice and their eggs on contact. The electric charge is so small that users will not notice it; it is absolutely safe to use. Comb activity is indicated by an audible signal.

The unit is powered by conventional AA batteries.

The **MEDISANA** electric lice comb **LCS** is suitable for adults and children over the age of 3. The comb can also be applied to household pets, such as dogs and cats.

### 3.1 Inserting / removing the battery

Remove the battery compartment lid **3** by lifting out the catch on the top end of the cover with your fingernail and removing it from the housing. Insert the included battery (type AA LR6). Observe the proper polarity (+/-) of the batteries as indicated on the batteries themselves and on the inside of the battery compartment. Replace the battery compartment lid **3**, by inserting the hinge on the lower end and then pressing the catch on the top end into the housing until it snaps in. Follow the same procedure when replacing a battery.

We recommend that you do not use rechargeable batteries. Remove the battery from the unit if you do not plan on using it again anytime in the near future.

#### BATTERY SAFETY INFORMATION

- Keep away from children!
- Do not short-circuit!
- Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop!
- Do not recharge!
- Do not throw into fire!

### 3.2 Operation

1. The electric lice comb **LCS** may only be used on dry hair. If necessary, detangle your hair and undo any knots or braids prior to using the comb.
2. To switch off the unit, slide the ON/OFF switch **1** forward/upward. A high-pitched humming sound indicates that the unit is in operation.
3. Position the comb at a 45° - 90° angle and guide it slowly through your hair.
4. While you are combing your hair, the hum will briefly stop whenever the teeth make contact with a louse or egg. If the humming sound does not return, the comb teeth are **4** clogged. Switch the comb off. Use the small cleaning brush **2**, that is integrated in the end of the handle, to clean the teeth **4**. Continue combing your hair. Allow enough time to thoroughly comb your hair.
5. Don't forget to switch off the unit after you are done and to clean it thoroughly. When not in use, always put the protective cap **5** back on the unit.
6. Since the **MEDISANA** electric lice comb **LCS** is particularly effective against adult lice, continue using the comb for about 2 weeks to ensure that you eliminate all lice.
7. To prevent the return of lice, we recommend that you use the electric lice comb once a week as a precautionary measure.

### 3.3 Helpful tips if you are afflicted by head lice

- You can supplement the application of your electric lice comb by using a chemical agent. For unusual cases, please consult a physician.
- Be aware that medicines are only available from a pharmacist and may only be available with a prescription. Due to their possible side effects, chemical agents should only be used if you have an actual affliction of head lice and not as a means of prevention. A physician can verify this and prescribe treatment.
- Head lice can survive on textiles for a brief period of time. For that reason, you should clean all objects that come into contact with your head. This measure can only be regarded as a supplement to treating affected persons. You cannot prevent the return of lice by "only washing or cleaning".
- If caps, scarves, bed linen, pillow cases, towels or stuffed animals come into contact with lice, they should be washed at 60 °C or higher or put into a dryer for 10 minutes. In you are unsure or are not able to wash these objects, you can enclose them in a hermetically sealed plastic bag for two weeks at room temperature or else place them into a freezer for one day at -15 °C.

### 4.1 Cleaning and maintenance

- Regular cleaning and good maintenance ensures optimum results and extends the life of your unit.
- Clean the lice comb carefully after each use, particularly if it is used by more than one person.
- Remove the batteries before you clean the unit.
- To clean the housing, use a dry, soft paper or fabric cloth. To clean the comb teeth, use the included cleaning brush **2**.
- Do not use any strong detergents for cleaning, e.g. scouring cream or other cleansing agents, as they may damage the surface of the unit.
- When cleaning, never immerse the unit in water and make sure that no water gets into the unit.
- Do not expose the unit to direct sunlight and protect it from dirt and moisture.
- Remove the batteries from the unit if you do not plan on using it again anytime in the near future. Otherwise, the batteries may start leaking.
- Store the unit in a clean, dry and cool location.

### 4.2 Disposal

This product must not be disposed together with the domestic waste.

All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.

Please remove the battery before disposing of the equipment. Do not dispose of old batteries with your household waste, but dispose of them at a battery collection station at a recycling site or in a shop.

Consult your municipal authority or your dealer for information about disposal.

### 4.3 Specifications

Name and model	: <b>MEDISANA</b> electric lice comb <b>LCS</b>
Dimensions (L x W x H)	: 173 x 63 x 23 mm
Weight	: approx. 75 g incl. battery
Power supply	: 1,5 V= , Battery (type AA, LR 6)
Article number	: 41000
EAN number	: 4015588410000

CE 0297

**In accordance with our policy of continual product improvement, we reserve the right to technical and design changes without notice.**

## 5.1 Warranty and repair terms

Please contact your dealer or the service centre in case of a claim under the warranty. If you have to send in the unit, please enclose a copy of your receipt and state what the defect is.

The following warranty terms apply:

1. The warranty period for **MEDISANA** products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
2. Defects in material or workmanship will be removed free of charge within the warranty period.
3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
4. The following is excluded under the warranty:
  - a. All damage which has arisen due to improper treatment, e.g. non-observance of the user instructions.
  - b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
  - c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
  - d. Accessories which are subject to normal wear and tear.
5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

### **MEDISANA AG**

Itterpark 7-9

40724 Hilden

Tel.: +49 (0) 2103 / 2007-60

Fax: +49 (0) 2103 / 2007-626

eMail: [info@medisana.de](mailto:info@medisana.de)

Internet: [www.medisana.de](http://www.medisana.de)

The service centre address is shown on the last page.

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D&L Marketing Oy
ALERE Electronics AB